PermaCem Universal

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Instructions for use English

Product description

PermaCem Universal is a universally applicable, dual-curing, bonding composite luting cement for the permanent cementation of indirect restorations and root posts. The radiopaque material is used together with a suitable bonding agent (e.g. LuxaBond Universal). The material is available in the Smartmix-Syringe.

Intended use

- Definitive luting of crowns, bridges, inlays and onlays made of metal (precious and non-precious metal), composite, PEEK or ceramic (oxide and silicate ceramics, for example zirconium oxide or lithium disilicate)
- Permanent luting of porcelain or composite veneers
- Permanent luting of crowns and bridges on implant abutments
- Cementation of root posts

Limitations of use

Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.

Indications

Dental cases requiring indirect restoration:

- Major tooth defects that cannot be restored by direct restoration, e.g. extensive carious lesions, traumatic injuries or extensive dental wear
- Endodontically treated teeth that are susceptible to fractures
- Tooth loss
- Dental misalignments and deformities, e.g. individual or multiple diastemata and internal discoloration
- It is not possible to predict the functional or esthetic outcomes of direct restorations
- Occlusal modification

Contraindications

- Do not apply PermaCem Universal directly to the exposed dental pulp
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended users

Dentist

Notes for use

- With the Smartmix version the material is applied directly without a dispensing device.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions for the light-curing unit.
- It is advisable to use a rubber dam.
- Bring the material to room temperature (15 – 25 °C/59 – 77 °F) before use.
- Ensure that the surface to be treated is not dirty before each following working step.
- During polymerization of the cement, an oxygen inhibition layer develops in the areas that have been exposed to atmospheric oxygen. This layer does not set. For this reason, only large areas of excess material should be removed, allowing smaller areas of excess material to remain. The oxygen inhibition layer is removed together with the smaller areas of excess material when finishing the edges.
- Observe manufacturer's information for other materials that are used together with PermaCem Universal.

Using the Smartmix-syringe (for image, see »Handling«)

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it (Fig. 1).
2. Attach a new mixing tip [B] (Fig. 2).

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix-syringe [C].

3. Stop the mixing tip by turning it 90° clockwise (Fig. 3).
4. The material is mixed in the tip when it is expressed and it can be applied directly.

Note: When using a syringe for the first time, squeeze out a pea-sized amount of material and discard it!

Note: After application, leave the used mixing tip on the Smartmix-Syringe as a cap!

Recommended use

Overview of the entire application

A.	Prepare tooth structure			
B.	Where necessary determine cement shade with Vitique Try In			
C.	Pretreat restoration materials			
	I. Pretreat oxide ceramic / metal:	II. Pretreat etchable silicate ceramics:	III. Pretreat Composite or PEEK restoration:	IV. Pretreat fiberglass-reinforced root posts (e.g. LuxaPost):
	Sand-blast, clean, apply primer	Etch, rinse with water, dry, apply Silane	Sand-blast, clean, apply primer (according to the manufacturer's instructions)	Rinse with water (according to the manufacturer's instructions)
D.	Apply cement to the pretreated tooth substance or restoration:			Apply cement into the pretreated root canal:
	Apply bond to the prepared tooth, apply cement, remove excess material, cure cement, if necessary light-cure cement			

A. Prepare tooth substance

1. Prepare the tooth substance in accordance with the general rules. Treat root canal endodontically and prepare for the root post insertion, leaving at least 4 mm guttapercha apically.
2. Remove any residual temporary luting material if necessary.
3. Rinse preparation with water and dry with oil-free air.
4. Try in the restoration.

B. Where necessary determine cement shade with Vitique Try In

Note: The cement shade is selected prior to placing the rubber dam.

1. Determine shade of the luting cement with Vitique Try In as per manufacturer's instructions.
2. If the chosen shade is corrected: Clean the prepared tooth and restoration with water and repeat Try In with a suitable shade.
3. Remove the restoration. Clean the restoration and prepared tooth thoroughly using a water spray.

Note: Any residual Vitique Try In Paste can affect the permanent cementation of the restoration.

Note: After contamination with Vitique Try In Paste clean the restoration with water spray and dry with oil-free and water-free air. Silanize silicate ceramic again if necessary (e.g. with Vitique Silane or LuxaPrimer).

C. Pretreat restoration materials

C.1. Pretreat oxide ceramic / metal

Note: Do not etch oxide ceramic.

1. Prepare the clean and dry restoration in accordance with manufacturer's instructions. If not otherwise described, sandblast the surfaces to be bonded with aluminum oxide (≤50 µm).
2. Clean the restoration with alcohol and dry with air that is free of oil and water.
3. Pretreat restoration with a suitable primer (e.g. LuxaPrimer or LuxaBond Universal) as per manufacturer's instructions.

C.II. Pretreat etchable silicate ceramics

1. Prepare the clean and dry restoration in accordance with manufacturer's instructions. If not otherwise described, rinse the surfaces to be bonded thoroughly with 5% liquid acid as per manufacturer's instructions.
2. Dry the restoration with oil-free and water-free air.
3. Prepare the surfaces to be bonded with a silanizing agent (e.g. LuxaPrimer, LuxaBond Universal or Vitique Silane) as per manufacturer's instructions.

C.III. Pretreat Composite or PEEK restoration

Note: Do not etch composite or PEEK restorations. Only use Isopropanol for cleaning PEEK.

1. Prepare the clean and dry restoration in accordance with manufacturer's instructions. If not described elsewhere, sand-blast the surfaces to be bonded with aluminum oxide until they have a matt appearance. Clean the restoration with alcohol and dry with air that is free of oil and water. Afterwards apply a primer (for composite e.g. LuxaPrimer or LuxaBond Universal; for PEEK Luxatemp Glaze & Bond).

C.IV. Pretreat fiberglass-reinforced root post (e.g. LuxaPost)

Pretreat in accordance with the manufacturer's instructions.

Note: Silanization is not necessary with the use of the LuxaPost.

D. Apply cement onto the pretreated tooth substance or restoration

Note: The use of a rubber dam to isolate the working area is recommended.

1. Apply a suitable adhesive (e.g. LuxaBond Universal) as per manufacturer's instructions.
2. Apply PermaCem Universal inside the pretreated restoration.

Note: For veneer cementation, it is recommended that a maximum of two veneers be cemented simultaneously due to the short curing time.

3. Insert restoration within 60 seconds after start of mixing and place onto the prepared tooth with light pressure and stabilize until the cement has cured.
4. Remove cement residues within 30 to 60 seconds after inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler.
5. Remove interdental excess carefully using dental floss.
6. Alternatively, the cement residues can be briefly (1–2 seconds) cured with light and removed.
7. Once mixing begins, allow the material to chemically cure completely for 07:00 min.
8. In the case of light-permeable restorations, additionally cure the material with a suitable light unit for at least 20 seconds.

D. Apply cement into the pretreated root canal

Note: The use of a rubber dam to isolate the working area is recommended.

Note: With use of PermaCem Universal in conjunction with LuxaBond Universal, the working time in the root canal is 30 seconds. The use of other bonding agents can reduce the working time.

1. Apply a suitable adhesive, ideally a dual-curing one, (e.g. LuxaBond Universal) as per manufacturer's instructions.
2. Apply PermaCem Universal directly into the prepared root canal with the aid of the corresponding tip. Introduce the tip as deeply into the root canal as possible.

Note: During the application of PermaCem Universal, ensure that the tip remains inserted in the material and is only moved upwards when the material moves upwards.

3. Insert the root post within 30 seconds after mixing has begun and cement using slight pressure.
4. Remove cement residues within 30 seconds after inserting the root post with the aid of a microbrush, brush, foam pellet or scaler.
5. Once mixing begins, allow the material to chemically cure completely for 07:00 min.

Note: The material can be additionally cured with a suitable light unit for at least 20 seconds.

Timing

Working time (including mixing time)	~ 60 seconds*
Working time in the root canal (including mixing time)	~ 30 seconds**
Removal of excess material (without light-curing)	30 - 60 seconds after the restoration is placed
Light-curing time for removal of excess (optional)	1-2 seconds
Maximum intraoral curing time (including mixing time)	07:00 min
Light-curing time for curing (optional)	Min. 20 seconds per side

*The time given applies to a room temperature of 23 °C. Higher temperatures can reduce this time, lower temperatures can increase this time.

**The working time in the root canal depends on the bonding agent used, as this can reduce the working time.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings / precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Performance characteristics of the device

Complies with ISO 4049:2019, Type 2, class 3 (incl. shade stability and radiopacity).

Storage / disposal

- Store in a dry place at 2 to 25°C (36 to 77°F)!
- After application, leave the used mixing tip on the syringe as a cap!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Acrylates and methacrylates (bis- and multi-functional, approx 29%): EBPADMA (approx. 9%), UDMA (approx. 8%), TEDMA (approx. 7%), HEMA (approx. 3%), barium glass/silicon dioxide (3 µm; approx. 69%), additives.

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