Icon

Inhalt | Contents

Gebrauchsinformation – DEUTSCH	4
Instructions for use – ENGLISH	14
Mode d'emploi – FRANÇAIS	23
Istruzioni d'uso – ITALIANO	34
Instruccion de uso – ESPAÑOL	44
Instruções de uso – PORTUGUÊS	54
Gebruiksaanwijzing – NEDERLANDS	63
Brugsanvisning – DANSK	73
Användarinstruktioner – SVENSKA	83
Instrukcja użycia – POLSKI	92
Информация по применению – РУССКИЙ	103
使用说明 – 中文	113
取扱説明書 – 日本語	122
사용 정보 – 한국어	130

Produktbeschreibung

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Product description

Icon is an innovative product for the micro-invasive treatment of cariogenic lesions in proximal regions and on smooth surfaces and available in two variations:

Icon Caries Infiltrant – Proximal is a specially developed hard tissue preserving treatment for incipient proximal caries.

Icon Caries Infiltrant – Smooth Surface is a specially developed infiltration product for smooth surfaces and particularly well-suited for orthodontic patients after bracket removal.

Indications

- Micro-invasive treatment of incipient caries (lesion depth up to D1) in proximal regions, see »Lesion depth classification*« on page 14.
- Micro-invasive treatment of smooth surface enamel lesions (e. g. cariogenic white spots after bracket removal).

Lesion depth classification*

* Radiographic lesion depth classification according to bitewing images.

Contraindications

Do not use this material:

- for deeper seated lesions (D2 – D3) (see »Lesion depth classification*« on page 14) or on cavitated enamel (enamel defects).
- in case of known allergies to any material component or existing contact allergies.

General safety information

- ▶ For dental use only!
- ▶ Keep out of reach of children!

Caution! Icon-Etch can cause chemical burns. Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.

- ▶ Limit the application to the surface to be treated. Avoid contact of the etching gel with the gingiva and adjacent teeth.
- ▶ Use only the enclosed syringe and application tip to apply the etching gel in the patient's mouth.
- ▶ Discard the application tip after use and cap the syringe.

Caution! Avoid contact of Icon-Infiltrant and Icon-Dry with eyes and skin.

- ▶ In case of eye contact, rinse with copious amounts of water, and, if necessary, consult a physician.
- ▶ In case of skin contact, wash the areas thoroughly with water and soap.

Caution! Improper use or incorrect application of Icon-Etch may cause pain in patients.

- ▶ Use Icon-Etch on enamel surfaces only. Do not apply onto exposed dentine or root cementum.

Please note

- Light-curing units should have an output of 450nm and should be checked regularly.
- The light intensity should be at least 800mW/cm².
- Place the light unit as close to the material as possible.
- On superficially remineralized, mostly older white spots the pseudo-intact surface layer may be fairly thick. In this case, a 2 minute etching step to remove the surface layer may not be sufficient and etching and rinsing may have to be repeated. The teeth are then dried with Icon-Dry and infiltrated.
- The etched surface must not be touched or contaminated with saliva or moisture until the treatment resumes. If contamination occurs after drying re-etch for approx. 10 seconds and repeat drying with Icon-Dry.
- Use a rubber dam. For best treatment results a sufficiently dry working area is critical. Therefore, placement of a rubber dam is imperative. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whale-dent/Hygienic). For smooth surface applications a liquid dam (Liquid Dam) can be used as well.

Please follow the manufacturer's specifications and instructions for use of the rubber dam.

- Discard syringes after use.

Proximal surface treatment

- The syringes included in the treatment unit contain sufficient material for treating two proximal lesions.
- Additional proximal surfaces can be treated consecutively during one visit.

Smooth surface treatment


- The syringes included in the treatment unit contain sufficient material for treating two to three smooth surface lesions.
- Multiple smooth surface lesions can be treated at the same time.
- For smooth surface applications a liquid dam (Liquid Dam) can be used as well.
- Treating orthodontic patients after bracket removal: Unless white spots are treated early, i. e. 1 to 2 months after bracket removal, it is recommended to repeat the etching step. A third etching process is recommended if a white spot is still visible after Icon-Dry is applied.
- In total, a lesion can be etched up to three times for 2 minutes each with Icon-Etch.

- For topical cariogenic white spots an area of 2 mm beyond the lesion site is etched. At the discretion of the dentist the entire smooth surface area ought to be etched and infiltrated in case of large area white spots as they occur after bracket removal.

Recommended use


Proximal surface treatment with Icon Caries Infiltrant – Proximal

1. Before the start of treatment, clean the affected tooth and adjacent teeth. Remove any cleaning residue with water spray. Place a rubber dam.


2. Introduce one of the enclosed dental wedges into the interdental.

- ▶ For better access to the proximal region the wedge handle can be bent or detached with a twisting motion. In order to obtain sufficient separation of the teeth the wedges must be seated at an adequate depth in the interdentalium.
- ▶ A separation of approx. 50 μm is recommended to ensure a successful treatment. This can be done gradually: Insert the wedge to the point of resistance. Wait 3 to 5 seconds and gradually move the wedge deeper until sufficient separation is achieved.
- ▶ Leave the wedge in the proximal space during the entire treatment procedure.


3. Screw the Proximal-Tip onto the Icon-Etch syringe, and introduce the application tip into the interdentalium. Be sure to align the green side of the Proximal-Tip with the area to be treated.


- ▶ **The material will dispense only on the green side of the application tip.**


4. Apply an ample amount of Icon-Etch onto the lesion site (1 ½ to 2 turns of the shaft). Let Icon-Etch set for 2 minutes. Remove excess material.
- ▶ **The screw syringe can be operated with one hand.**


5. Remove the application tip from the interdentium. Suction off Icon-Etch and rinse with water for a minimum of 30 s. Dry with oil-free and water-free air.

6. Screw the application tip onto the Icon-Dry syringe, apply an ample amount of material onto the lesion and let set for 30 s. Dry with oil-free and water-free air.

- ▶ **Do not apply Icon-Infiltrant under direct operating light as this may cause the material to set.**
- 7. Screw a new Proximal-Tip onto the Icon-Infiltrant syringe, and introduce the application tip into the interdentium. Be sure to align the green side of the Proximal-Tip with the area to be treated.
- ▶ **The material will dispense only on the green side of the Proximal-Tip.**
- 8. Apply an ample amount of Icon-Infiltrant onto the lesion site (1 ½ to 2 turns of the shaft).
- 9. Let Icon-Infiltrant set for 3 minutes and add material if necessary.
- 10. Remove the application tip from the interdentium. Remove excess material with dental floss.


11. Light-cure Icon-Infiltrant from all sides for at least 40 seconds (total).

12. Screw a new Proximal-Tip onto the Icon-Infiltrant syringe, repeat the application (steps 7 – 8), and let set for 1 minute. Remove the application tip from the interdental. Remove excess material with dental floss. Then light-cure from all sides for at least 40 seconds (total).
13. If a second lesion needs to be treated repeat steps 1 through 12 for the affected site.
14. Remove the wedge and rubber dam. Use polishing strips for the surface finish.


Smooth surface treatment with Icon Caries Infiltrant – Smooth Surface

1. Before the start of treatment, clean the affected tooth and adjacent teeth. Remove any cleaning residue with water spray. Place a rubber dam.
2. Screw the Smooth Surface-Tip onto the Icon-Etch syringe.


3. Apply an ample amount of Icon-Etch onto the lesion site by turning the syringe shaft carefully, and let set for 2 minutes. Remove excess material with a cotton roll.
- **The screw syringe can be operated with one hand.**

Note: Unless white spots are treated early, i. e. 1 to 2 months after bracket removal, it is recommended to repeat the etching step. A third etching process is recommended if a white spot is still visible after Icon-Dry is applied.


4. Suction off Icon-Etch and rinse with water for at least 30 seconds. Dry with oil-free and water-free air.


5. Screw the application tip onto the Icon-Dry syringe, apply an ample amount of material onto the lesion, and let set for 30 seconds. Dry with oil-free and water-free air.

Note: When wetted with Icon-Dry, the whitish-opaque coloration on the etched enamel should diminish. If this is not the case repeat the etching step once or twice for 2 minutes each, and rinse and dry the teeth again (steps 3-5).

6. Screw a new Smooth Surface-Tip onto the Icon-Infiltrant syringe.
 - ▶ **Do not apply Icon-Infiltrant under direct operating light as this may cause the material to set.**

7. Apply an ample amount of Icon-Infiltrant onto the etched surface by turning the shaft.
8. Let Icon-Infiltrant set for 3 minutes. Remove excess material with a cotton roll and dental floss.

9. Light-cure Icon-Infiltrant for 40 seconds.
10. Screw a new Smooth Surface-Tip onto the Icon-Infiltrant syringe, repeat the application (steps 6 – 7), and let set for 1 minute. Remove excess material with a cotton roll and dental floss, and light-cure for a minimum of 40 seconds.
11. Remove the rubber dam. Use polishing cups (or similar) for the surface finish.

Side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon-Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Interaction

Polymerization inhibiting substances such as products containing eugenol (e. g. temporary cements) must not be used in combination with this material.

Composition

- Icon-Etch: Hydrochloric acid, pyrogenic silicic acid, surface-active substances
- Icon-Dry: 99% ethanol
- Icon-Infiltrant: Methacrylate-based resin matrix, initiators, additives

Storage

- ▶ Store at 2 to 25 °C/35,6 to 77 °F!
- ▶ Protect from direct sunlight!
- ▶ Do not use after the expiration date!

Patient card + sticker

This material is not radiopaque. To document the infiltrated surfaces and lesion depths the enclosed patient card can be filled out at the time of treatment and follow-up visits and handed to the patient. In addition, the enclosed sticker supports the documentation in the patient's record.

1. Front: Enter the name of the patient and office stamp.
2. Tooth diagram: Mark the treated area.
3. Record the treated tooth [1], the treated surface [2] (d = distal, m = mesial, v = vestibular (buccal, cervical, labial), l = lingual or palatal), lesion depth*, and the treatment date [3].

1 Zahn | Tooth 2 Fläche | Surface

3 Läsionstiefe | Lesion depth

Behandlungsdatum | Treatment date

Kontrolldatum 1 | Follow-up appointment 1

Kontrolldatum 2 | Follow-up appointment 2

Kontrolldatum 3 | Follow-up appointment 3

E1 E2 D1 D2 D3 P

4. Follow-up visit [4]: Record the lesion depth* and date of the follow-up appointment. Annual check-ups are recommended.

see »Lesion depth classification« on page 14.

Packaging

Icon Caries Infiltrant - Proximal	
Starter Kit @ 2 treatment units	REF 220341
Package @ 7 treatment units	REF 220237

Icon Caries Infiltrant – Smooth Surface	
Starter Kit @ 2 treatment units	REF 220343
Package @ 7 treatment units	REF 220238

One treatment unit includes:

Icon Caries Infiltrant - Proximal

- 1 Icon-Etch syringe @ 0.3 ml fluid
- 1 Icon-Dry syringe @ 0.45 ml fluid
- 1 Icon-Infiltrant syringe @ 0.45 ml fluid
- Accessories

Icon Caries Infiltrant – Smooth Surface

- 1 Icon-Etch syringe @ 0.45 ml fluid
- 1 Icon-Dry syringe @ 0.45 ml fluid
- 1 Icon-Infiltrant syringe @ 0.45 ml fluid
- Accessories

Caution: Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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The diagram illustrates the patient card and marking sheet. It consists of two main parts: 'Zahn | Tooth' and 'Fläche | Surface'. The 'Zahn' part shows a tooth diagram with a vertical line indicating the treatment area. The 'Fläche' part shows a square diagram with a vertical line indicating the treatment area. The 'Zahn' part has a table for recording treatment data, including 'Läsionsstiefe | Lesion depth', 'Behandlungsdatum | Treatment date', and 'Kontrolldatum 1 | Follow-up appointment 1'. The 'Fläche' part has a table for recording treatment data, including 'Läsionsstiefe | Lesion depth', 'Behandlungsdatum | Treatment date', and 'Kontrolldatum 1 | Follow-up appointment 1'. The card also includes a section for 'Kontrolldatum 2 | Follow-up appointment 2' and 'Kontrolldatum 3 | Follow-up appointment 3'.

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The diagram illustrates the components of the patient pass and the dental chart. On the left, a tooth is shown with a lesion depth indicated by a vertical line. The tooth is labeled 'Zahn | Tooth' and the surface is labeled 'Fläche | Surface'. The patient pass is shown on the right, with fields for 'Lesionsdjup | Lesion depth', 'Behandlingsdatum | Treatment date', and 'Kontrolldatum 1 | Follow up appointment 1', 'Kontrolldatum 2 | Follow up appointment 2', and 'Kontrolldatum 3 | Follow up appointment 3'. The dental chart shows a tooth with a lesion depth and a treatment date. The chart is divided into sections for 'E1', 'E2', 'D1', 'D2', 'D3', and 'P'.

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The diagram illustrates the components of the patient booklet and the corresponding sticker. The booklet includes a section for 'Läsionstiefe | Lesion depth' with a grid for recording lesion depth (E1, E2, D1, D2, D3, P) and a section for 'Behandlungsdatum | Treatment date'. The sticker shows a tooth diagram with markings for 'Zahn | Tooth' and 'Fläche | Surface'. It also includes a timeline for follow-up appointments (Kontrolldatum 1, 2, 3) with a magnifying glass icon indicating the need for control.

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DMG Chemisch-Pharmazeutische Fabrik GmbH