OsteoBio® Apatos - Italiano -

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OsteoBio® Apatos - English -

OsteoBio® bone matrix, Class III Medical Device (ref. Dir. 93/42/EEC amended by 2007/47 EEC) for bone grafting. CHARACTERISTICS: carbonated nanocrystal bone mineral and organic matrix of natural heterologous origin, obtained by tissues treated by a patented exclusive TecnoSS process. Osteoconductive; resorbable; biocompatible; it can function as a carrier for selected medication and drugs. Available in: a) granules of cancellous bone, cortical bone and as proper mix of both; b)pre-shaped block of natural cancellous bone.

GRANULOMETRY: 0.6 – 1 mm (Cortical bone/Cancellous bone/Mix); 1 - 2 mm (Mix). Product is slightly radiopaque.

CLINICAL INDICATIONS: collagenated filler of non infected, non sclerotic, non subjected to loading and blood bedewed bone defects

PACKAGING: glass vial or disposable single use blister.

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use. In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: OsteoBio® Apatos in granules must be hydrated before grafting with sterile saline and/or with patient blood and/or with the selected drug for surgery. Moreover, to improve handling, OsteoBio® Apatos in granules can be mixed with OsteoBio® Gel 0: the mix obtained can be grafted in the receiving site through a sterile spatula or biomaterial syringe in sub-periosteum site. Grafting site must be well-vascularised and the product must be properly stabilized and protected. OsteoBio® Apatos Sp-Block must be hydrated before grafting with sterile saline or with antibiotics; then it could be adapted to grafting site, which has to be carefully curretled to guarantee a proper blood bedewing of grafting. Product must be stabilized in grafting site through osteosynthesis microscrews and it is recommended to be protected from risk of exposure through a resorbable membrane/barrier. Grafting site must be well-vascularised and the product must be properly stabilized and protected.

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and +40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): AXXXXFY, ASXXXXFY, ACXXXXFY, AXy, ASXy, ACXy, AXyXX, ASXyXX, ACXyXX, ABLXXy. The letter "Y" indicates the origin of material (S= porcine; E= equine).

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