OsteoBio[®] Derma - Italiano –

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OsteoBio[®] Derma- English –

OsteoBio[®] membrane: Class III Medical Device (ref. Dir. 93/42 CEE amended by Dir.2007/47 EEC) for tissue regeneration.

CHARACTERISTICS: collagenated membrane from porcine derma, obtained with a patented and exclusive Tecnos process. Its consistency and resistance allow a perfect stabilization and long lasting protection of graft before if flaps are properly and carefully sutured. Resorbable; biocompatible; it can function as a carrier for selected medication and drugs.

SIZES: Standard, 2-3 mm thick ; Fine, 1 mm thick ; X-Fine 0.4 - 0.6 mm.

ESTIMATED RESORPTION TIMES: Standard about 4/6 months; Fine: about 2/4 months; X-Fine: about 1 months. Resorption of material may require longer time, depending on grafting site characteristics.

CLINICAL INDICATIONS: Membrane for the protection of reconstructions or partial or total recoveries of lost tissue portions (laparocoele); collagenated filler for bone distraction and split-crest procedures; protection of non infected, non sclerotic and well blood bedewed bone defects

PACKAGING: single use sterile blister.

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use. In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: before use membrane must be hydrated by dipping it in lukewarm sterile saline for 3-5 minutes. Once product has reached the desired plasticity, it can be shaped with sterile instruments and placed in the grafting site properly prepared. In Odontostomatology, it is recommended to prepare a pocket with an elevator in order to achieve a perfect hold of the membrane in site after stitching the flaps. NB: In the event of exposure, membrane should be removed from grafting site only in case of supra-infection , because the consistency and plasticity of product allow a second intention healing of wound.

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and +40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): EMXXZY. The letter "Y" indicates the origin of material (S= porcine).

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