OsteoBio® Evolution - Italiano –

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OsteoBio® Evolution - English -

OsteoBio® membrane: Class III Medical Device (ref. Dir. 93/42/CEE amended by Dir. 2007/47 EEC) for tissue regeneration.

CHARACTERISTICS: high tightness and resistance; high stability but soft and flexible consistency;over 2 months resorption times; easy hydration; possibility to function as a carrier of medication and selected drugs; possibility to suture or stick with fibrin glue; tenacious and resistant to traction; it provides a valid space-making able to guarantee a regeneration of underlying tissues, without infiltration and infection risks.

SIZES: Standard (about 0.6 mm thick). Fine version (about 0.4 mm thick). X-Fine (about 0.2 mm thick). Evolution Perio kit 6 pre-shaped periodontal membranes blister (0.3 mm thick).

ESTIMATED RESORPTION TIMES: Standard about 4/6 months; Fine about 3/4 months; X-Fine and Perio-kit about 2 months. Resorption of material may require longer time, depending on grafting site characteristics.

CLINICAL INDICATIONS: Evolution membrane is suitable for all surgical procedures requiring a membrane able to guide a perfectly vital and vascularised connective tissue neo-formation and/or filling or replacement of connective tissues (General Surgery, Thoracic Surgery, Neurosurgery, Gynaecology; Maxillofacial Surgery, Otorhinolaryngology, Plastic Surgery, Reconstruction Odontostomatology Surgery).

PACKAGING: single use sterile blister.

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use. In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: membrane can be manipulated with sterile instruments, properly cut and shaped and then eventually hydrated with antibiotics or other selected drugs. Once it reaches the desired shape and plasticity, it must be placed in the receiving site and stabilized through suture thread or osteosynthesis fixtures or pins or by sticking it through fibrin glue in order to better take advantage of its resistance and stability properties.NB: In the event of exposure, membrane should be removed from grafting site only in case of supra-infection , because the consistency and plasticity of product allow a second intention healing of wound.

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and + 40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): EVXXLZE, EVXXHZE, EVOLZE, EVOHZE, EMXXZY,EMXXZS,EM9XZS. The letter "Y" indicates the origin of material (S= porcine; E= equine).

PRESEVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and + 40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): EVXXLZE, EVXXHZE, EVOLZE, EVOHZE, EMXXZY,EMXXZS,EM9XZS. The letter "Y" indicates the origin of material (S= porcine; E= equine).

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