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OsteoBio® Lamina | English

OsteoBio® bone matrix, Class III Medical Device (ref. Dir. 93/42/CEE amended by 2007/47 EEC) for bone grafting.

Lamina: cortical bone matrix in pre-shaped bone lamina.

CHARACTERISTICS: carbonated nanocrystal bone mineral and collagen of natural heterologous origin, obtained by cortical bone tissue treated by a patented exclusive Tecnos process. Natural substitute of autogenous or allogous bone, it conserves the same intimate structures (matrix, porous form); osteoconductive; resorbable; biocompatible; hygroscopic; it can function as a carrier for selected medication and drugs.

SHAPES (to be verified with distributors): NORM: medium mineralization grade; stiffness of pre-shaped similar to natural bone. SOFT: decreased mineralization grade; flexible pre-shaped.

CLINICAL INDICATIONS: co-adjvant for the reconstruction or the partial or complete recovery of lost bone portions, as fillers of non infected, non sclerotic and well blood bedewed bone defects.

PACKAGING: single use sterile blister.

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use.

In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: before use, Lamina must be hydrated by dipping it in lukewarm sterile saline for 3-5 minutes. Once product has reached the desired plasticity, it can be placed in the grafting site properly prepared. For the codes LS10HSL/LS10HE do not hydrate before use. Lamina can be adapted to the anatomical characteristics of receiving site by sterile scissors or other instruments and then it must be stabilized in grafting site through osteosynthesis microscrews or directly through suture to surrounding tissues with a non-traumatic triangular section needle. In the event of exposure, Lamina should be removed from grafting site only in case of supra-infection , because the consistency and plasticity of product allow a second intention healing of wound. Grafting site must be well-vascularised and the product must be properly stabilized and protected. Product must be placed in sub-periosteum site.

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and +40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): LNXXYZ, LSXXZY. The letter "Y" indicates the origin of material (S=porcine; E=equine).

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