#### OsteoBio® mp3 - Italiano –

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CONTROINDICAZIONI: non utilizzare in pazienti che presentano patologie o affezioni contraindicating surgery o presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

GRANULOMETRY: 0.6 – 1 mm; 2.4 mm ( hospital marketing only). Product is slightly radiopaque.

CLINICAL INDICATIONS: collagenated filler of non infected, non sclerotic, non subjected to loading and well blood bedewed bone defects.

PACKAGING: disposable single use sterile syringes

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use. In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: product is ready-to-use and no hydration is needed before grafting. It can be directly grafted in receiving site. After grafting, receiving site must be closed by flap repositioning and suture. Product can function as a carrier of selected medication and drugs. In Oral Surgery, mp3® must be grafted in sub-periosteum site.

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and +40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): A3XXXFYA2XXXFY. The letter "Y" indicates the origin of material (S=porcine; E=equine).

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