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OsteoBio® Putty - English –

OsteoBio® bone matrix, Class III Medical Device (ref. Dir. 93/42/CEE amended by 2007/47 EEC) for bone grafting. CHARACTERISTICS: paste of bone and collagen of natural heterologous origin. The collagen component is represented by OsteoBio® Gel 0, type I and III collagen matrix with Omega-3 and oleic-linoleic poly-unsaturated fat acids, with anti-oxidant effect on free radicals and with a valid aid in tissue regeneration. Heterologous bone component provides the scaffold and regenerative properties, thus making the product suitable for bone and tissue regeneration interventions. Moreover, product can be used as a carrier of hydro soluble and/or liposoluble drugs.

GRANULOMETRY: ≤ 0,300 mm; product is slightly radiopaque.

CLINICAL INDICATIONS: Odontostomatology: alveolar filler for immediate post-extractive implants, where it enhances primary stability; treatment of dehiscences and peri-implantitis; bone filler after tooth extraction and granulomas, dentigenous cysts and split crest procedures. Orthopaedy and Traumatology: filler of cages for spinal surgery; carrier of granules bone fillers and drugs.

PACKAGING: disposable single use sterile syringes

COLLATERAL EFFECTS: none known so far. Although the manufacturing process proved the inhibition of adverse antibody response, nevertheless this adverse response cannot be totally excluded.

CONTRAINDICATIONS: use is advised against in patients with pathologies or affections contraindicating surgery or presenting acute or chronic infections (osteomyelitis) in progress. Use is not recommended in patients affected by immunologic disorders or auto-immune pathologies, in particular children or elderly. Use is not recommended directly or in proximity to big blood vessels or nerves. Use is not recommended in scarcely blood bedewed, necrotic or infected grafting sites. Use is not recommended in sites not allowing the correct stabilization and protection of product once grafted.

WARNING: the product has not been experimented on pregnant patients so far.

Product designed for surgeons with solid experience in biomaterial grafting: please strictly follow the guidelines supplied by the scientific literature on xenografts, in addition to these instructions for use. In each application, the specialized surgeon is responsible for the proper clinical management of biomaterials and the patient's care.

HANDLING: product is ready to use and can be taken out with a sterile spatula and placed directly in receiving site. Due to its plastic consistency, OsteoBio® Putty must be grafted only in cavities able to contain it in stable way, thus it is not recommended in one or two-wall defects or in sinus lifting with lateral access procedures. If necessary, it is possible to mix the product with granules bone fillers or with drugs. In Odontostomatology it is recommended to protect the grafting with a resorbable membrane, in order to increase primary stability, before proceeding with suture of integumental flaps. Grafting site must be well-vascularised, product must be placed in sub-periosteum site after proper stabilization and protection of grafting. N.B.: to fill big cavities (vol > 1 cm3), one part of product should be mixed with three parts of OsteoBio® Gen-Os or OsteoBio® Chips (volume ratio of Putty: 25%).

EXPIRY DATE: the validity of product is 5 years from sterilization date and expiry date is impressed both in labels and on box. Product must not be used if packaging is open or damaged; the product is sterile and disposable: re-sterilisation is absolutely prohibited.

PRESERVATION: product must be preserved from direct contact with hot surfaces and solar rays, and stored in dry environment between +4 and + 40°C temperature.

ORIGIN OF PRODUCT: in labels and on boxes, one of the following codes is impressed (REF): HPTXXY. The letter "Y" indicates the origin of material (S=porcine; E= equine.

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