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P.G.A. Aragó

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User Instructions

P.G.A. Aragó

Braided and Coated Polyglycolic Acid
Surgical Suture, sterile and absorbable

Description:
PGA. Aragó is a sterile, synthetic, braided and absorbable surgical suture based on polyglycolic acid (PGA). The empirical formula is (C₄H₆O₂). In order to enhance the sliding and handling properties, the braided thread has been covered by a mixture composed of stearic acid with a glycolic and lactic acid copolymer (30/70) or sacroic fatty acid esters and polycaprolactone. For better visibility, the thread has been dyed with D&C Violet No. 2 dye (C.160725), but it can be supplied colourless. It complies with the requirements of the European Pharmacopoeia for Sterile Absorbable Braided Synthetic Sutures and the American Pharmacopoeia (USP) with respect to Absorbable Synthetic Sutures.

Indications:
PGA. Aragó is indicated for applications where, in accordance with surgical practice, absorbable materials are used for sutures. It is suitable for the approximation of soft tissue and ligatures, but not for cardiovascular or neurological tissue.

Effect:
PGA. Aragó causes a minimal acute reaction of the tissue on which a gradual encapsulation occurs due to the connective tissue. The gradual loss of resistance to breakage and eventual absorption of the PGA. Aragó is carried out through hydrolysis. The absorption begins with a loss of resistance to tension, followed by a loss of mass. Studies on implants in rats show the following profile:

After two weeks, the remaining resistance is higher than 80 % of the initial resistance. At three weeks, the remaining resistance is approximately 50-55 % of the initial resistance depending on the USP number. Total dissolution is completed at approximately 90 days.

Contraindications:
Due to the fact that the PGA. Aragó is absorbable, it must not be used on injuries that depend on prolonged reinforcement with sutures.

Warnings:
PGA. Aragó users must be familiar with surgical techniques and the use of absorbable sutures, because the risk of the injury reopening (dehiscence) depends on where it is used and the suture material employed. Circumstantially, the PGA. Aragó is not suitable for elderly or weak patients, or for those suffering from a slow healing of wounds. Since this is a resorbable suture, the surgeon must consider using a supplementary non-absorbable suture in the suture of parts that may be subject to expansion, tension or distension or require additional support. Do not re-sterilize/re-use. Do not use if the packaging is open. Re-using a single-use product creates a potential risk of infection for the patient. Product contamination may cause injury or disease, or lead to the death of the patient. Re-processing could jeopardise the structural integrity of the product or cause it to fail. Using the product without the label and instructions will result in using the product beyond its expiry date and the loss of its traceability. Altering the original packaging can jeopardise the sterility and functionality of the product.

Precautions:
Sutures that must remain on the skin for more than 7 days may cause local irritations, for which they must be cut or removed. In tissues with deficient blood supply, there could be a rejection of the thread caused by the slow absorption. All sutures must be handled with

care. Avoid damaging the suture by crushing or folding it due to through the use of forceps or needle holders. If the needles are remodelled manually, they may lose their strength and become less resistant to bending and breaking. The surgeon must take into consideration the type of knot to be used in a surgical technique. When handling sutures, due precautions must be taken to prevent surgical instruments, such as forceps and needle holders, from damaging the suture material by crushing, folding or scraping it, and/or deforming the needle. The user must be familiar with suture techniques and must handle the needles with care so as to not be pricked by them. Discard used needles into a special container. Do not use a suture between one-third (1/3) and a half (1/2) of the distance from the distal extremity to the point of the needle. If it is held close to the tip, the capacity for penetration will decrease and this could also cause a fracture to the needle. If the needle is held by the distal extremity, the needle may bend or fracture.

Side effects:
Minor reactions in the tissue surrounding the suture material may appear at the beginning. As with any foreign body, the PGA. Aragó may exacerbate a prior infection. As with any foreign body, prolonged contact with this or any other suture with saline solutions such as those found in urinary and biliary conduits could lead to the formation of calculi.

Presentations:
This suture can be supplied in USP sizes 10/0 (EP 0.2) up to USP 6 (EP 8) in pre-cut threads of different lengths and with a variable number of threads, without a needle or swaged to different types of atraumatic needles. See the catalogue for more detailed information.

Storage conditions:
Store at room temperature, prevent prolonged exposure to high temperatures and protect from dampness. Do not use after the expiry date.

Explanation of symbols on the packaging:

	Batch No
	Catalogue reference
	Single use only
	Do not re-sterilise
	Expiry yyyy-mm-year/month
	See instructions
	Sterilised with ethylene oxide
	Keep out of sunlight
	Do not use if the unitary packaging is open or deteriorated
	Manufacturer

CE Marking + identification number of the notified body. This product meets the essential requirements of Council Directive 93/42/EEC, concerning medical devices.

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