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Poliglecaprona Aragó

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User instructions

Poliglecaprone Aragó

Mono-filament glycolic and caprolactone Surgical Suture, sterile and resorbable

Description:
Poliglecaprone Aragó is a sterile, synthetic, mono-filament and resorbable surgical suture/ligation which is a glycolic and ε-caprolactone copolymer with low tissue reactivity. The empirical molecular formula is C₁₄H₂₀N₂O₇. For better visibility, the strand has been coloured using the D&C Violet N° 2 (C.I.60725) although it is available colour-free. Complies with the monographs of European Pharmacopoeia and American USP for mono-filament resorbable synthetic sutures, except for small variations in diameter.

Indications:
Poliglecaprone Aragó is indicated for applications where, in accordance with surgical practice, absorbable materials are used for sutures. It is suitable for the approximation of soft tissue and ligatures, but not for cardiovascular or neurological tissue.

Effect:
Poliglecaprone Aragó causes a minimal acute reaction of the tissue to which a gradual encapsulation occurs due to the connective tissue. The gradual loss of resistance to breakage and eventual absorption of the Poliglecaprone Aragó is carried out through hydrolysis. The absorption begins with a loss of resistance to tension, followed by a loss of mass. Studies on implants in rats show the following profile:

After 1 week, the residual resistance is approximately 60%-75% of the initial resistance. After 2 weeks, the residual resistance is approximately 40% of the initial resistance depending on the USP number. Complete dissolution takes approximately 90-110 days.

Contraindications:
Due to the fact that the Poliglecaprone Aragó is absorbable, it must not be used on injuries that depend on prolonged reinforcement with sutures.

Warnings:
Poliglecaprone Aragó users must be familiar with surgical techniques and the use of absorbable sutures, because the risk of the injury resuming (dehiscence) depends on where it is used and the suture material employed. Circumstantially, the Poliglecaprone Aragó is not suitable for elderly or weak patients, or for those suffering from a slow healing of wounds. Since this is an absorbable suture, the surgeon must consider using a supplementary non-absorbable suture in the suture of parts that may be subject to expansion, tension or distension or requiring additional support. Do not re-sterilize/re-use. Do not use if the packaging is open. Re-using a single-use product creates a potential risk or infection for the patient. Product contamination may cause injury or disease, or lead to the death of the patient. Re-processing could jeopardise the structural integrity of the product or cause it to fail. Using the product without the label and instructions will result in using the product beyond its expiry date and the loss of its traceability. Altering the original packaging can jeopardise the sterility and functionality of the product.

Precautions:
Sutures that must remain on the skin for more than 7 days may cause local irritations, for which they must be cut or removed. In tissues with deficient blood supply, there could be a rejection of the thread caused by the slow absorption. All sutures must be handled with care.

Avoid damaging the suture by crushing or folding it through the use of forceps or needle holders. If the needles are remodelled manually, they may lose their strength and become less resistant to bending and breaking. The surgeon must take into consideration the type of knot to be used in a surgical technique. When handling sutures, good precautions must be taken to prevent surgical instruments, such as forceps and needle holders, from damaging the suture material by crushing, folding or scraping it, and/or deforming the needle. The user must be familiar with surgical techniques and must handle the needles with care so as to not be pricked by them. Discard used needles into a special container. Hold the needle between one-third (1/3) and one-half (1/2) of the distance from the distal extremity to the tip of the needle. If it is held close to the tip, the capacity for penetration will decrease and this could also cause a fracture to the needle. If the needle is held by the distal extremity, the needle may bend or fracture.

Side effects:
Mild reactions in the tissue surrounding the suture material may appear at the beginning. As with any foreign body, Poliglecaprone Aragó may exacerbate a prior infection. As with any foreign body, prolonged contact with this or any other suture with saline solutions such as those found in urinary and biliary conduits may lead to the formation of calculi.

Presentations:
This suture can be supplied in sizes USP 7/0 (EP 0.5) up to USP 2 (EP 5), in pre-cut threads of different lengths and with a variable number of threads, without a needle or swaged to different types of traumatic needles. See the catalogue for more detailed information.

Storage conditions:
Store at room temperature, prevent prolonged exposure to high temperatures and protect from dampness. Do not use after the expiry date.

Explanation of symbols on the packaging:

	Batch No
	Catalogue reference
	Single use only
	Do not re-sterilise
	Expiry yyyy-mm-year/month
	See instructions
	Sterilised with Ethylene Oxide
	Keep out of sunlight
	Do not use if the unit packaging is open or deteriorated
	Manufacturer

CE Marking = identification number of the notified body. This product meets the essential requirements of the Council Directive 93/42/EEC, concerning medical devices

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