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### Seda Trenzada Aragón

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#### User Instructions

### Braided Silk Aragón

#### Multifilament braided silk Surgical, sterile and non-absorbable sutures

**Description:** Braided Silk Aragón is a sterile, braided and non-resorbable surgical suture comprised of braided strands of silk fibres (polypeptide chain), obtained from Bombyx mori L silkworm cocoons. For better visibility, the strand has been dyed black with dye Sulpho Black 1 (C.I. 53185) or Logwood Black (C.I. 75290). It is also available without dye (white). The suture material is covered with a mixture of wax and/or silicone. Meets the requirements for non-resorbable sterile fibres according to European Pharmacopoeia and American USP.

**Indications:** Braided Silk Aragón is indicated for application where non-resorbable suture materials are used for surgical procedures. It is suitable for closing soft tissues and for ligations, including use for general, skin, oral and ophthalmology surgery.

**Effect:** Braided Silk Aragón causes a minimal acute reaction in tissue with gradual encapsulation due to the connective tissue. It slides through tissue easily. It is not resorbable. Provides temporary support to the wound until the healing reaches a degree allowing it to resist normal tension. Although it is classified as "non-resorbable material", it is digested through proteolysis after an extended period of time.

**Contraindications:** Due to the loss of resistance which can occur after long periods of "in vivo" implantation, silk sutures are not recommended when resistance must be maintained permanently.

**Warnings:** Braided Silk Aragón users must be familiar with surgical techniques and the use of non-absorbable sutures, because the risk of the injury reopening (dehiscence) depend on where it is used and the suture material employed. Do not re-sterilise/re-use. Do not use if the packaging is open. Re-using a single-use product creates a potential risk of infection for the patient. Product contamination may jeopardise integrity or lead to the death of the patient. Re-processing could jeopardise the structural integrity of the product or cause to fail. Using the product without the label and instructions will result in using the product beyond its expiry date and the loss of its traceability. Altering the original packaging may jeopardise the sterility and functionality of the product.

**Precautions:** All sutures must be handled with care. Avoid damaging the suture by crushing or folding it through the use of forceps or needle holders. If needles are remodelled manually, they may lose their strength and become less resistant to bending and breaking. The surgeon must take into consideration the type of knot to be used in a surgical technique. When handling sutures, due precautions must be taken to

prevent surgical instruments, such as forceps and needle holders, from damaging the suture material by crushing, folding or scraping it and/or deforming the needle. The user must be familiar with suture techniques and must handle the needles with care so as to not be pricked by them. Discard used needles into a special container. Hold the needle between one-third (1/3) and a half (1/2) of the distance from the distal extremity to the tip of the needle. If it is held close to the tip, the capacity for penetration will decrease and this could also cause a fracture to the needle. If the needle is held by the distal extremity, the needle may bend or fracture.

**Side effects:** Mild reactions in the tissue surrounding the suture material may appear at the beginning. As with any foreign object, Braided Silk Aragón may exacerbate a prior infection. As with any foreign object, prolonged contact of this or any other suture with saline solutions such as those in urinary and biliary conduits may lead to the formation of calculus.

**Presentations:** This suture can be supplied in USP sizes 10/0 (EP 0.2) up to USP 6 (EP 8) in pre-cut threads of different lengths and with a variable number of threads, without a needle or swaged to different types of traumatic needles. See the catalogue to obtain more detailed information.

**Storage conditions:** Store at room temperature, prevent prolonged exposure to high temperatures and protect from dampness. Do not use after the expiry date.

#### Explanation of symbols on the packaging:

	Batch No
	Catalogue reference
	Single use only
	Do not re-sterilise
	Expiry yyyy-mm- year/month
	See instructions
	Sterilised by Gamma Rays
	Keep out of sunlight
	Do not use if the unit packaging is open or deteriorated
	Manufacturer

**CE 0318** CE Marking + identification number of the notified body. This product meets the essential requirements of the Council Directive 93/42/EEC concerning medical devices

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