# PREDCLEAR

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# PREDCLEAR

## High transparent Matrix Material

### **EN** **Instruction for use** [MD] UE Medical device

#### **Product description**

PREDCLEAR is a high transparent matrix material based on addition curing vinyl silicones. For this purpose make an impression of the surface before beginning the preparation or after the model has been waxed up. This impression builds a tailor made transparent matrix for the subsequent restoration.

#### **Indications/Intended use**

Matrix material especially for

- the composite injection technique
- temporary crown and bridge work
- the composite layering technique
- for the transfer of brackets or fibers from the model to the mouth

#### **Performance features**

The performance features of the product meet the requirements of the intended use.

#### **Contraindications**

None known.

#### **Patient target group**

Persons who are treated during a dental procedure.

#### **Intended users**

This medical device should only be used by a professionally trained dental practitioner.

#### **Preparing of the Automix cartridge**

Remove and dispose the closure-cap and replace it with a supplied mixing cannula. Mount the cartridge into the application gun and extrude needed material.

#### **Note:**

Store used cartridge with fixed used mixing cannula. Do not re-use the cartridge cap. Closing the cartridge by turning on the original cap between use could cause cross contamination of components and thus premature setting.

#### **Application**

Load a non-perforated impression tray with the automatically mixed PREDCLEAR and insert the tray not later than 45 seconds from start of mixing (working time at 23 °C or 74 °F) into the mouth.

Maintain position until the material has completely set (about 2:00 minutes) and remove the matrix from the mouth.

#### **Disinfection**

Disinfection of the matrix with disinfection solutions specific for materials based on addition curing silicones is possible. Refer to the manufacturer instructions.

#### **Storage**

Do not store above 25 °C (77 °F). Protect from direct

sunlight.

Do not use after expiry date.

#### **Additional notes:**

- Do not use latex gloves because they may inhibit the polymerization of vinyl silicones. Do not touch tooth with latex gloves.
- Keep away from children!

#### **Disposal**

Disposal of the product according to local authority regulations.

#### **Reporting obligation**

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to Laboratoires PRED (pharma@pred.fr), the manufacturer and the competent authority.

#### **Warranty**

First Scientific Dental Materials GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental Materials GmbH's sole obligation shall be repair or replacement of the First Scientific Dental Materials GmbH product.

#### **Limitation of Liability**

Except where prohibited by law, First Scientific Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.