

ANA 2000[®]

High Copper Admixed Alloy

CAPSULES

NEW CAPSULES



ANA 2000[®]

High Copper Admixed Alloy

Our disposable capsules have been optimized reaching yet a higher level of quality. The new special design gives a homogenous mix. Being self activating and now possible to open without a rod they are as convenient as possible. With a precise dosage in each disposable capsule the consistency is guaranteed – consistently!

Today when time is of the essence it is important to have a dispensing system that allows quick and easy handling, the ANA 2000 capsules with low mercury exposure makes this possible.



ANA 2000 CAPSULES

- No. 1 – Green capsule
(400 mg alloy + 400 mg mercury)
- No. 2 – Blue capsule
(600 mg alloy + 600 mg mercury)
- No. 3 – Grey capsule
(800 mg alloy + 800 mg mercury)

CAPSULES

■ FACT	● BENEFIT
New design	Optimal amalgamation
Pre-dosed	Exact consistency
	Time saving
Self activating	Easy to mix
Colour coded	Easily differentiated



The capsules are available in three different sizes and are delivered in boxes of 50 or 450 capsules.



COMPOSITION

Ag 43,1%
Sn 30,8%
Cu 26,1%

PHYSICAL PROPERTIES

Compressive strength:
1 h 140 MPa
24 h 450 MPa
7 days 510 MPa

Static creep: 0.1%

Dimensional change:
24 h +0.02%

ANA 2000

■ FACT	● BENEFIT
High copper	Long-lasting, resistant and strong
Controlled expansion	Excellent marginal adaptation
Lathe cut particles	Crunchy feeling when condensing
Short setting time	Time saving



Box 1082, SE-262 21 Ängelholm, Sweden
Int tel: (+46) 431 44 33 60 - Fax: (+46) 431 44 33 99
mail@nordiskadental.se - www.nordiskadental.se
Certified acc. to ISO 9001 and ISO 13485