

**Tokuyama**

# **Universal Bond II**

Self-cured Dental Universal Adhesive

The data contained in this report are obtained in Tsukuba Research Laboratory in Tokuyama Dental Corporation. The data are measurement representative values, not guaranteed values.

# Tokuyama

# Universal Bond II

## Self-cured Dental Universal Adhesive



# Difference from TOKUYAMA UNIVERSAL BOND

## 1. Storage temperature

- Can be stored at room temperature  
*(previous UB version needed to be refrigerated)*
- Storage at a temperatures between 0 to 25° C (32 to 77° F )

## 2. Color changed

- Can be visually confirmed that it is mixed by color of the liquid



Tokuyama  
**Universal Bond II**  
Self-cured Dental Universal Adhesive

# Difference from TOKUYAMA UNIVERSAL BOND

## 3. Volume of Bond B reduced

- The volume of Bond B of Tokuyama Universal Bond II is 4mL (5mL for Tokuyama Universal Bond)



The volume of one drop is adjusted as "Bond A: Bond B = 1: 0.8".

# Product Information

## Tokuyama Universal Bond II Kit




1 × Bond A (5mL)



1 × Bond B (4mL)



1 × mixing well



5 × 1-Dimple Disposable mixing wells



25 × Disposable applicator brushes