

Alveolar Ridge Preservation with **Benex[®] II**

What does it mean?


Alveolar Ridge Preservation means the treatment of the dental alveolus after extraction. 3 months after the **Benex[®]** extraction you find a considerably better ridge relation than with conventional gentle extractions.

The **Benex[®]** finds its successful application in both, private practices and universities. The **Benex[®]** has achieved an excellent status worldwide as basis for a subsequently successful implantation.

The new support for the dismantled **Benex[®]** System in a washbasket guarantees an optimal cleanability of **Benex[®]** in a washing machine or in an ultrasonic bath. All components can be fixed safely in the support.

Upon cleaning, the system can be sterilized in assembled condition.

In addition to the **Benex[®]** components, an optional periotome, an Xtool and an optional **Benex[®]** Pole Extractor can be placed in the support.


12.302.00


Benex[®] II Extraction System,
consisting of:
see ♦ in the table of
page 05-04/05

12.303.00


Benex[®] II Basic Kit,
consisting of:
see ♦ in the table of
page 05-04/05



Optional Instrumentary


12.300.11

Benex[®] Pole Extractor


12.300.45

Benex[®] Blade

47.525.51

Driver Guide FD

47.525.55

Ratchet with demountable handle

Illustration	Article Description	Order Quantity
	12.302.00 Benex [®] II Extraction System, consisting of: Benex [®] II Extractor, 2 Pullropes 48 mm, Driver Guide, Screw short 1.6 mm + 2.1 mm, Screw long 1.6 mm + 2.1 mm, 1 Drill ea. for 1.6 mm, 2.1 mm Screws, Quadrant Support, 85.195.10 Washbasket with Lid, 12.302.01 Tray / Rack for Benex [®] II	1 piece
	12.302.01 Benex [®] II Tray / Rack	1 piece
	85.195.00 Washbasket 1/1 with Lid and Press Button Lock	1 piece
	12.303.00 Benex [®] II Basic Kit, consisting of: Extractor, 2 Pullropes 48 mm, Driver Guide, Screw short 1.6 mm + 2.1 mm, Screw long 1.6 mm + 2.1 mm, 1 Drill ea. for 1.6 mm + 2.1 mm Screws, Quadrant Support, 85.194.10 Washbasket with Lid	1 piece
	85.194.10 Washbasket 1/2 with Lid and Press Button Lock	1 piece
	12.300.08 Benex [®] II Extractor	1 piece
	12.300.15 Replacement Support Disc, 8 mm (PTFE), optional	1 piece
	12.300.16 Support Disc, diagonally right, optional	1 piece
	12.300.17 Support Disc, diagonally left, optional	1 piece
	12.300.20 Pullrope, 48 mm	2 pieces

Illustration	Article Description	Order Quantity
	12.300.30 Diamond coated Drill for Screws Ø 1.6 mm 12.300.60 and 12.300.70	2 pieces (1  )
	12.300.35 Diamond coated Drill for Screws Ø 2.1 mm 12.300.65 and 12.300.75	2 pieces (1  )
	12.300.47 Driver Guide, short	1 piece
	12.300.60 Screw, Ø 1.6 mm, 10 mm, S = Short	2 pieces (1  )
	12.300.65 Screw, Ø 2.1 mm, 10 mm, SF = Short & Fat	2 pieces (1  )
	12.300.70 Screw, Ø 1.6 mm, 16 mm, L = Long	2 pieces (1  )
	12.300.75 Screw, Ø 2.1 mm, 16 mm, LF = Long & Fat	2 pieces (1  )
	12.300.80 Quadrant Support for Benex ® to bridgeover bigger gaps and for the universal molding	1 piece
	12.300.11 Benex ® Pole Extractor, optional	1 piece
	47.525.55 Ratchet with demountable handle, optional	1 piece
	47.525.51 Driver Guide FD, optional	1 piece
	12.300.45 Benex ® Blade for Driver Guide FD, optional	1 piece

